



DO IN WATER

Oxygen is poorly soluble in water. The amount of oxygen in water depends on physical, chemical and biochemical activities taking place in water body. The solubility of DO in water at saturation at any temperature & pressure is given by Henry's law. The saturation value of DO in water is of the order of 8 to 15 mg/l depending upon the temperature & pressure. At 20°C temperature and standard pressure the maximum amount of oxygen that can dissolve in fresh water is 9 ppm. If the temperature decreases, there may be more oxygen dissolved in the sample water.

PRACTICAL RELEVANCE

DO level	Water quality
0-4.0 ppm	Poor, undesirable odours, Effect aquatic life
4.1-7.9 ppm	Fair
8.0-12 ppm	Good
12+	Retest

1. Fall in DO levels causes undesirable odours, tastes and reduce the acceptability of water for domestic uses.
2. In aquatic life, as DO drops below 4 ppm fish and other species are affected. When DO level is below 2 ppm fishes are threatened and may get decline.
3. Determination of DO is important in industrial purposes. In steam generation, DO is important factor causing corrosion control of the boiler materials.



DETERMINATION OF DISSOLVED OXYGEN BY IODOMETRIC TITRATION METHOD

I. SAMPLING

FOR DO determination-

- (i) The sample container should be made of Glass/ BOD bottle. The sample water should be taken in a stoppered bottle very carefully without any air bubbles which could raise oxygen level by aerating the sample.
- (ii) Minimum sample size should be 250ml and type Grab.
- (iii) The sample should be analyzed immediately as recommended within 15 minutes.

Location:


Date & Time of Sampling

Sam p. No.	Sample Described as	Sampling Source	Container	Sample Size & Type	Preservation	Max. Storage Recommended


II. THEORY

Dissolved oxygen (DO) is determined by (Winkler's method) IODOMETRIC TITRATION method. However since dissolved oxygen in water is in molecular state and is not capable of reacting with *KI*, therefore an oxygen carrier is to be produced, which enables the dissolved oxygen to take part in reaction.

The method introducing adding of Conc. Soln. of $MnSO_4$ and alkaline *KI* into the water sample. A white ppt of $Mn(OH)_2$ which is oxidized by dissolved oxygen present in water to give a brown ppt. of $MnO(OH)_2$ which acts as an oxygen carrier to enable the dissolved oxygen to take part in reaction.



This $MnO(OH)_2$ in acidic medium dissolves to produce $[O]$ which is also in acidic medium Oxidizes *KI* and liberates free I_2 in an equivalent amount of dissolved oxygen present in the water sample.



This liberated I_2 is then titrated against standard hypo solution using starch as an indicator.



METHOD	
---------------	--

III. PREPARATION

Apparatus: stoppered bottle 250 ml, pipette 5 ml, conical flask, pipette 50ml, burette 50 ml, syringe 5ml.

Chemicals: $MnSO_4$, Alkaline *KI*, Conc. H_2SO_4 , $Na_2S_2O_3$, $5 H_2O (N/40)$, starch indicator.


- (i) $MnSO_4$ soln.: Dissolve 5 gm of $MnSO_4$ in d.w and make up the volume 100ml.
- (ii) Alkaline *KI* Soln.: Dissolve 40 g *NaOH* and 20 g *KI* in d.w and make up the volume to 100ml.
- (iii) Conc. H_2SO_4 Soln.: Take nearly about 50ml of Conc. H_2SO_4 in a bottle.
- (iv) $Na_2S_2O_3, 5 H_2O (N/40)$ soln.: Dissolve 6.20 g of $Na_2S_2O_3$ in d.w and make up the volume to 1lit.
- (v) Starch soln.: Dissolve 1 g powdered starch with d.w to form a paste and add 100 ml of boiled water to it with constant stirring, then boil for 5 minutes and then cool.




III. PROCEDURE

(A) Standardisation of $\text{Na}_2\text{S}_2\text{O}_3$ soln.

1. Rinse & fill up the burette with $\text{Na}_2\text{S}_2\text{O}_3$ soln.
2. Rinse & pipette out 20 ml std $\text{K}_2\text{Cr}_2\text{O}_7$ Soln and taken in a conical flask.
3. 3 ml KI soln (by Pipette) and then 3ml conc. H_2SO_4 soln (by syringe) are added to it and covered the mouth with watch glass. The mixture is Shaken well and now the soln become DEEP BROWN Colour (due to liberation of Iodine)




4. It is then titrated (liberated I_2) against $\text{Na}_2\text{S}_2\text{O}_3$ Soln until the brown colour fades to pale yellow colour.
5. Now, 3ml Starch Soln is added (by syringe), then the soln turns DEEP BLUE Colour. Continued the titration until blue colour disappears leaving a COLOURLESS Soln (End point)




(B) Determination of DO


1. Rinse & fill up the burette with stand $\text{Na}_2\text{S}_2\text{O}_3$ Soln.
2. A Known amount of Sample water (Say 250ml) is taken in a stoppered bottle.
3. 3 ml MnSO_4 Soln (by pipette) is added, dipping the end well below the surface of water. Also, 3ml Alkaline - Iodide-acid Soln (by pipette) is added and shaken well.
4. Appearance of brown ppt (of basic manganic oxide) indicates the presence of DO. The ppt is allowed to form and settle down.



5. 3ml conc H_2SO_4 Soln (by Syringe) is added to it and stoppered the bottle. The mixture is Shaken well to dissolve the ppt completely and now the soln become DEEP BROWN Colour (due to liberation of Iodine)



6. Rinse & measured out 100 ml of the above soln and taken in a conical flask.
7. It is then titrated (liberated I_2) against Stand $\text{Na}_2\text{S}_2\text{O}_3$ Soln until deep brown colour fades to pale yellow colour.
8. Now, 3ml starch indicator soln (by syringe) is added, when the Soln turns DEEP BLUE Colour. Continued the titration until blue colour disappears leaving a colourless Soln (End point)



Sl.No.	Vol.of Sample soln. Taken, ml V_1	Vol.of titrant $Na_2S_2O_3$ used, ml V_2			str. of $Na_2S_2O_3$ N_2	str. of water N_1
		initial	final	diff		
					$(N/40)$?

$$V_1 N_1 = V_2 N_2$$

(water) (hypo)

$$\Rightarrow N_1 = \frac{V_2}{50} \times \frac{1}{40} (N)$$

$$\begin{aligned} \text{Strength of dissolved Oxygen} &= \frac{V_2}{50 \times 40} (N) \times 8 \text{ g/l} && [\text{eq. wt. of } O_2 = 8] \\ &= \frac{V_2}{50 \times 40} \times 8 \text{ g/l} \times 1000 \text{ mg/L} \\ &= \text{mg/L} \end{aligned}$$

V. RESULTS & COMMENTS

Standard Value: as per	Experimental Value
Comments	





COLOURLESS

WATER SAMPLE



HAZY BROWN PPT.

BASIC MANGANIC OXIDE



CLEAR BROWN SOLN.

IODINE SOLN

$MnSO_4$
 $NaOH-KI-NaNO_3$

conc. H_2SO_4



PALE YELLOW SOLN

IODINE SOLN (low conc)

$NO_2S_2O_3$



Starch Indicator

BLUE-BLACK SOLN

STARCH-IODIDE SOLN



COLOURLESS SOLN

END-POINT

$Na_2S_2O_3$

DETERMINATION OF THE DISSOLVED O_2 IN WATER BY IODOMETRIC TITRATION METHOD

